

# Instructions for use (IFU)

Rev 1 / April 2024

Download the latest IFU from [www.hicareindia.com](http://www.hicareindia.com)

GTIN	Product	Description
8908019116304	Sterilised Camera Cover	Sterilised Laparoscopic Camera Cover 100"x 8.5" (2540 x 216mm)

## Indications for use

Hi-Care Remedy's, Sterilised Laparoscopic Camera Cover / Camera Sleeve is intended to provide a sterile interface between the operation theatre camera and the healthcare professional. The Camera Cover should be used by trained healthcare professionals.

The Camera Cover is sterile, for single patient/procedure use and disposable.

## Instructions for use

1. Pull the Camera Cover so that it covers the cable. Ensure that the lens of the Camera is not obstructed. Secure the Camera Cover using the provided self-adhesive sticker and elastic bands.
2. Inspect the Camera Cover to ensure there are no holes or tears.
3. After use, dispose the product in accordance with your facility's internal guidelines

## Storage Conditions

Store in a cool dry place.



Manufacturer



Batch code/ Lot number



Date of Manufacture



Use by date



Do not re-use/ Single use only



Consult instructions for use or consult electronic instructions for use



Medical Device



Latex Free



Do not use if package is damaged



Sterilised using ethylene oxide



Keep dry



Keep away from sunlight